

MEDICAL OEM EQUIPMENT	MODEL NUMBER / DESCRIPTION	LUXTEL CERALUX LAMP TYPE	LUXTEL LAMP MODULE PART NO.
Excelitas/PerkinElmer/ILC	Cermax 175W, PE175BF, ILC LX175, ORC 175F	CL175BF	CL1339
Excelitas/PerkinElmer/ILC	Cermax 300W, PE300BF, ILC LX300, ORC 300F	CL300BF	CL1338
Excelitas/PerkinElmer/ILC	ME300BF	CL1938	
Acufex	FRX300 Xenon, Ultra Lite	CL300BF	
Aculux	AX 3001, AX 3002	CL300BF	
Acuity	A09501/A09502 Ultralife#013845	CL300BF	
ACMI/Circon	ALV-1, MV-9086 & MV-9087	CL300BF	CL1338
ACMI/Gyrus ACMI	MV-9090, XLS-300	CL1583	
Applied Fiber Optics	Britelite III, #01-30-40	CL300BF	
Applied Fiber Optics	Sunbeam Xenon Light Source	CL1579	
Athrex	SCB300	CL300BF	
Baxter/V. Mueller	VS5500, VS5600, LA8000	CL300BF	
Baxter/Edwards	OPTX300	CL300BF	
BEI/Xylog	300W Xenon	CL300BF	
BFW	Maxenon Xi300, MID 3010	CL300BF	CL1428
Birtcher/Solos/Woodbine	ELS2, GS9205, GS9250, GS9270, GS9480, GS9481	CL175BF	CL1339
Birtcher/Solos	GS9272, GS9255	CL300BF	CL1338
CirCon/AMCI	ALV-1 & MV-9086	CL300BF	CL1338
Concept/Linvatec	8430, 8180, 8170	CL300BF	
Cooper Surgical	Endomed LSS-700	CL300BF	
CUDA Fiberoptics	M2-300, XLS 300	CL300BF	CL1338
Dyonics(Smith & Nephew)	Dyo-Bright 3000XE, 300XL	CL300BF	1407
Fujinon/Hoya	EPX-4450, EPX-4440, X-402, X-4400, EPS301, EPX301, EPX302, EPX304, EPX401, FLX2000, FLX300, FLX300A, FLX600	CL300BF	
Fujinon/Hoya	EPK-i	CL1576	
Gulf Fiberoptics	Visualux LX 300	CL300BF	CL1338
Kay	Swallowing Workstation	CL300BF	
Intuitive	da Vinci Surgical System	CL300BF-10F or CL1582	
Isolux	XSB1300 300W Illuminator, Millenium	CL300BF	CL1338
Isolux	XSB1180 180W Illuminator	CL175BF	CL1339
Leica	M525	CL1938	
Leica	M720F60		CL1585
Leica	OH5, OH4, OH3, OHS1, MS3, MS2, MS1, F40	CL300BF	CL1426
Leisegang	LMS-3000 Illuminator	CL300BF	
Leitz	Endoscope	CL300BF	
Lorenz	Lumenon Xenon Headlight	CL300BF	CL1338
Linvatec/ConMed	8430, 8180, 8170	CL300BF	CL1426
Linvatec/ConMed	LS7700	CL1579	CL1752
Luxtec/Integra	LX300, 6000, 7000	CL300BF	
Luxtec/Integra	9300XDP, 9300XSP	CL300BF	1739
Luxtec/Integra	9300	CL300BF	CL1401 (Lamp Drawer)
Luxtec/Integra	9175	CL175BF	CL1402 (Lamp Drawer)
Medical Dynamics	6100XE, 6600XE	CL300BF	
Moeller-Wedel	FS2-23		CL1681
Nextec	300W Xenon, 1300	CL300BF	
Olympus	XLS(175W) Y-1176 or CLV-A(175W)	CL175BF	
Olympus	CLVS, CV-1, CLV-10, CLV-F10, CLV-20, CLV-S20, CLV-U20, CLV-S30, CLV-S40, CLV-U40, CLV-U60, CLV-160, CLV-180, CLV-260, CLV-260SL, ILX6300	CL300BF	
Olympus CLV-190/CLV-290	CLV-190/CLV-290 (MAJ-1817)	CL1854	
Olympus CLV-180	CLV-180 (MD-631)	CL1846	
Origin Medical	300 Xenon	CL300BF	
ORC	ORC6000, 5800 & 5800E	CL300BF	CL1338
Pentax	EPK-i5000/5010 EPK-i7000/7010	CL1952	
Pentax	EPK-i	CL1576	
Pentax	EPM3500, EMP3000, 3300, LX3000	CL300BF	
Pilling Weck	300W Xenon, 521393	CL300BF	CL1338
Products For Medicine	White Sun 300W, XL300	CL300BF	CL1338
QED	3010 300W Light Source	CL300BF	CL1338
Smith & Nephew	500XL	CL1938	CL1710
Smith & Nephew	300 XL	CL300BF	CL1336
Snowden Pencer	Xenon L/S 300W	CL300BF	
Storz (175 Watt Systems)	615 B or C, 20131520, 20132120	CL175BF	
Storz	300(20133020) 20133206 Lamp Only	CL300BF	
Storz (175 Watt Module)	20132025	CL175BF	
Storz (175 Watt Lamp)	20132026	CL175BF	
Storz (300 Watt Systems)	487, 610, 611, 630, NOVA 20134001, D-Light C 20133620, 20133120, 20133020, 20134020	CL300BF	
Storz (300 Watt Modules)	20133625, 20133027	CL300BF	
Storz (300 Watt Lamp)	20133028	CL300BF	
Stryker X8000	Module Swap - send in old module with iButton		CL1950
Stryker	X6000	CL300BF	CL1287
Stryker	X7000	CL1492	CL1571 (no housing)
Stryker	X8000	CL1495	CL1567 (no housing)
Sunoptics Surgical	Titan 300		CL1736
Sunoptics Surgical	SolarMaxx 300	CL300BF	CL1338
Sunoptics Surgical	SolarMaxx 180, Titan 180	CL175BF	CL1339
Toshiba	TRE-3000	CL300BF	
Video Dynamics	XLC-300	CL300BF	
Weck (Pilling)	FRX 2000, WECKLITE 300, 521393	CL300BF	
Welch-Allyn	49820	CL300BF	CL1338
Welch-Allyn	Pro Xenon		CL1681
Wolf	5052	CL1137	CL1177
Wolf	300W Xenon, 5131, 5151XE	CL300BF	CL1338
Zeiss	Superlux 300, Superlux 301 (with timer)	CL300BF	CL1426
Zeiss	Superlux 175 (with timer)	CL175BF	CL1429

\*Subject to change without notice. Updated: April 10, 2018.