

Stirrup Comparison Chart

PRODUCT GUIDE	Kyra P/N Equivalent	Allen® P/N Equivalent	Patient Weight Rating	Larger Boot for 327kg?	Flexible Adjustable Fin?	Amount of Pressure Relieving Padding	Lithotomy Range of Motion?	VELCRO® Foot Closure System?	Easy to Clean Silicone Closure Straps?	Seamless Waterproof Boot Pad?	Completely VELCRO® Free?
Allen® Yellofin® Stirrups	KYRA3000	O-YFASI	350 lbs	NA	NO: rigid hard plastic boot that is not flexible	1.3 cm	+84 -33: 117 degrees	YES: VELCRO® foot closure	NO: uses VELCRO®	NO: boot pad has seams	NO: pad attached to boot with VELCRO®
Kyra Comfort™ 350 Stirrup with traditional Straps	KYRA3000	O-YFASI	350 lbs	NA	YES: flexible adjustable fin	2.0 cm	+85 -55: 140 degrees	YES: VELCRO® foot closure	NO: uses VELCRO®	NO: boot pad has seams	NO: pad attached to boot with VELCRO®
Kyra Comfort™ 350 Stirrup with Silicone Foot Closure Straps	KYRA3100	No real equivalent but closest to OYFASI	350 lbs	NA	YES: flexible adjustable fin	2.0 cm	+85 -55: 140 degrees	NO: uses cleanable heavy duty silicone straps	YES: easy to clean silicone straps	NO: boot pad has seams	NO: pad attached to boot with VELCRO®
Allen® Yellofin® Elite	KYRA5000	O-YFES	500 lbs	NO: the 350 and 500 lb Yellofins® have the same size boot	NO: rigid hard plastic boot that is not flexible	1.3 cm	+84 -33: 117 degrees	YES: VELCRO® foot closure	NO: uses VELCRO®	NO: boot pad has seams	NO: pad attached to boot with VELCRO®
Kyra Comfort™ 500 Stirrup with traditional Closure straps	KYRA5000	O-YFES	500 lbs	YES: larger flexible boot sized to 95th percentile 327 kg patient	YES: flexible adjustable fin	2.5 cm	+85 -55: 140 degrees	YES: VELCRO® foot closure	NO: uses VELCRO®	NO: boot pad has seams	NO: pad attached to boot with VELCRO®
Kyra Comfort™ 500 Stirrup with Silicone Closure straps	KYRA5100	No real equivalent but closest to OYFES	500 lbs	YES: larger flexible boot sized to 95th percentile 327 kg patient	YES: flexible adjustable fin	2.5 cm	+85 -55: 140 degrees	NO: uses cleanable heavy duty silicone straps	YES: easy to clean silicone straps	NO: boot pad has seams	NO: pad attached to boot with VELCRO®
Kyra Clean™ 500 Stirrup with Silicone Closure Straps and Easy to Clean Seamless Boot Pad	KYRA5200	No equivalent	500 lbs	YES: larger flexible boot sized to 95th percentile 327 kg patient	YES: flexible adjustable fin	2.5 cm	+85 -55: 140 degrees	NO: uses cleanable heavy duty silicone straps	YES: easy to clean silicone straps	Seamless waterproof boot pad	YES: uses patent pending attachment system



build something better

888-611-KYRA (5972)

kyramedical.com



T +613 9550 1819

M + 61 409 865 105

E Info@Endovision.com.au

W Endovision.com.au