



ABC-ULTRA

FROM THE MANUFACTURER OF PINKIES



ABSORBABLE BED COVER - ULTRA

Made from a non-woven PAA material, with a quilt of high fluid absorbing SAP crystals and a blue PE backing, to prevent fluids escaping from beneath.

With an absorbable width of 75cm, (20cm more than most other commercially available pads) and the ability to absorb up to 12 litres of fluids, the ABC Ultra is superior in performance in all critical fluid management environments.

FEATURES

Waterproof backed blue PE material, kind to the skin so the patient feels dry

The super absorbent "quilt" absorbs up to 12 Litres of fluids

Single use, disposable and hygienic, saves on laundry bills

APPLICATIONS

Patient under pad (full bed coverage)

Operating Theatre table cover for extensive fluid absorption eg. Post Partum Haemorrhage & Caesarean Section

Suitable for use with bariatric patients

Protection of areas vulnerable to damaging secretions

Spill sheet and worktop work area covers

Suitable for Operating Theatre, Maternity, Emergency, Wards etc.

ABC-ULTRA - 225cm x 103cm - Box 20 \$110 + GST*

*Price at 09/2020

Our focus is patient care



ADVOCATES
FOR PATIENT
SAFETY AND
QUALITY
EQUIPMENT